

MIT OpenCourseWare
<http://ocw.mit.edu>

7.13 Experimental Microbial Genetics

Fall 2008

For information about citing these materials or our Terms of Use, visit: <http://ocw.mit.edu/terms>.

CAS Agar for Pseudomonas

Note: There is no Mg(II) in the medium

1. Medium, in a 1 l flask

540 ml H₂O

1.62 g sodium succinate

1.68 g casamino acids

0.43 g Na₂SO₄

5.44 g PIPES

Dissolve, adjust to pH=6.7 (exactly!) with 5 M NaOH

9 g agar

Autoclave (with the magnetic stirrer bar still in the flask)

2. CAS solution, in a 100 ml flask

30 ml H₂O

35 mg chrome azurol S

Dissolve, while stirring add

6.2 ml 1 mM FeCl₃ (made in 10 mM HCl)

3. Cetrимide-CAS solution

24 ml H₂O

47 mg cetrимide

Dissolve, then add the CAS solution slowly while stirring.

Autoclave

4. CAS plates

Cool medium to 50°C

Add Cetrимide-CAS solution slowly while stirring.

Pour for low iron CAS plates.

Add 1/1000 vol. of 100 mM FeCl₃ (made in 10 mM HCl) to a final conc. of 100 μM for high iron CAS plates.

UO/2/2001

Expected result for *Pseudomonas aeruginosa* strain PAO1 and *fur* mutants on high and low CAS medium plates

